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INSTRUCTIONS FOR USE

General Information & Intended Use

DURA-VIT Dental Implants are implantable medical devices used in case of total or partial mandibular and/or maxillary edentulism and/or dental agenesis. They get in contact with bone and mucous tissue to replace one or more missing natural teeth, up to the reconstruction of the entire dental arch. They are a support for a fixed prosthesis or an anchoring of removable prosthesis. The expected benefit is both functional, to rehabilitate correct and complete chewing, as well as to restore the ability to speak, and aesthetic, such as smiling.

The external design of the implants has a conical shape with a thread that allows to penetrate the harder bone, maintaining its stability in the softer bone, without damaging the most sensitive areas of the mouth as the mandibular canal and the Schneider's membrane (maxillary sinus). Therefore, the threads of the implants are different depending on the indications of use and the biological structures, that the patients present. Implants are provided of internal connection that allows the devices intended for prosthetic rehabilitation, to be anchored to them for the restoration of the missing elements.

The connection of the lines DURA-VIT 3P, EV, WIDE and PTERYGO feature an internal hexagon, platform switching and a Morse taper connection at 5° in total. The connection is unique for all the lines: any component, such as the abutment, can be placed into the implants regardless of implant diameter.

While the connection of the line SLIM feature an internal hexagon and a 22° conical connection, therefore the components to be used differ from the other lines.

The external surface of the implants is rough for promote the osseointegration, which is obtained with a double etching treatment for the lines DURA-VIT 3P, EV, WIDE and SLIM Ø3.4, or the sandblasting for the lines SLIM Ø3 and PTERYGO.

Material in use

Dental implants of the lines DURA-VIT 3P, EV, WIDE and SLIM Ø3.4 are realized in pure medical titanium of grade 4, with a treated surface with double etching.

Dental implants of the lines DURA-VIT SLIM Ø3.4 and PTERYGO, as well as the cover screws inside of the primary packaging, are supplied in titanium alloy Ti6Al4V (grade 5).

Just the transmucosal screws are produced in titanium alloy Ti6Al4V (grade 5) and peek.

Materials results compliant with harmonized standards.

In particular, the allergy to titanium is a very rare but possible event, so it is always necessary to check in advance with patients, who do not have this type of allergy.

Supply

DURA-VIT dental implants are available in sterile primary packaging.

The primary packaging of 3P, EV, WIDE and PTERYGO lines is a blister containing: 1 implant, 1 cover screw and 1 transmucosal screw.

The primary packaging of DURA-VIT SLIM is a blister containing: 1 implant, 1 cover screw.

The blister is inserted into a protective cardboard box (secondary packaging).

Each case contains one blister, traceability labels, package leaflet and implant card.

Traceability labels report the product code or product description, batch number and the UDI. They are stuck to the blister, Tyvek side, and they can be easily removed and applied to the implant card given to the patient and to the patient file for future reference.

Available diameters and lengths by line

DURA-VIT Implant line	Diameter (mm)	Length (mm)						
		Ø	6.5	8.0	10.0	12.0	14.0	16.0
SLIM	3.0	N/A	✓	✓	✓	✓	N/A	N/A
	3.4	N/A	✓	✓	✓	✓	N/A	N/A
3P	3.75	N/A	✓	✓	✓	✓	N/A	N/A
	4.0	✓	✓	✓	✓	✓	N/A	N/A

	4.5	✓	✓	✓	✓	✓	N/A	N/A
	5.0	✓	✓	✓	✓	✓	N/A	N/A
EVOLUTION	4.0	✓	✓	✓	✓	✓	✓	N/A
	4.5	✓	✓	✓	✓	✓	N/A	N/A
	5.0	✓	✓	✓	✓	✓	N/A	N/A
WIDE	5.5	✓	✓	✓	✓	✓	N/A	N/A
	6.0	✓	✓	✓	✓	✓	N/A	N/A
	6.5	✓	✓	✓	✓	✓	N/A	N/A
	7.0	✓	✓	✓	✓	✓	N/A	N/A
PTERYGO	4.7	N/A	N/A	N/A	N/A	N/A	✓	✓

*N/A: Not available

Warning

- Under law, only surgeons and certified dentists are allowed to use the medical devices in surgery. Any use by other persons is strictly prohibited. Therefore, DURA-VIT dental implant system must be handled and placed only by qualified and specially trained personnel.
- B. & B. Dental assumes no responsibility for any loss or damage/injury to the patient or any serious health hazards caused through improper use, or any use by unqualified or inadequately trained personnel, that is not in accordance with this instruction manual.
- The manufacturer excludes any responsibility for any damage to the patient from the use of medical devices not manufactured and not supplied by B. & B. Dental. Itself.
- Medical devices covered by these IFU must be placed in a hospital or clinical environment.
- Please check "Table 1" for a better understanding of the symbols on medical device labels.
- Please check "Table 2" for a better understanding of the symbols on the implant card.
- Potential adverse events associated with the use of this product may include, but are not limited to, the following:
 - Fractured implant
 - Rounding of the internal hexagon
 - Failure of or inadequate osseointegration
 - Unintended perforation of bone and/or soft tissue, such as inferior alveolar canal, maxillary sinus, and lingual artery, major palatine foramen.
 - Fracture of labial or lingual plate
 - Partial nerve injury or irreversible damage (paraesthesia/hyperaesthesia/anaesthesia)
 - Loss of bone due to concomitant events, such as local or systemic infection, peri-implantitis, loss of soft tissue volume

Proper selection of patient candidates, implant placement and post-operative care, combined with correct product use may help to mitigate these risks.

General contraindications

Dental implants and related accessories shall not be placed in case of **non-compliant patient** and if there is any **intraoral contraindication**, or if there is clinical or radiographic evidence of local pathology. Chronic or acute infectious diseases, chronic maxillary sinus osteitis, vascular problems, systemic diseases shall also be evaluated.

Absolute contraindications

Known allergies to titanium, recent myocardial infarction, valvular prosthesis, severe kidney disorder, severe liver disorders, osteomalacia, insulin-resistant diabetes, recent high-dosage radiotherapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.

Relative contraindications

Chemotherapy, moderate kidney disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged corticosteroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Precautionary information

- Since the surgical techniques are highly specialised and complex procedures, it is strongly recommended that practitioners shall complete a training course to learn established techniques of placement in oral implantology. Failure to seek proper training or follow instructions for use (IFU) may impact clinical success of implant, such as failure of the implant, loss of bone, and postoperative complications, that could lead to pain for the patient and seriously endanger his or her health.
- Surgical placement and consequent restoration of dental implants require intricate procedures that may involve risks. It is the responsibility of the practitioner to inform the patient of any risk(s) and alternative treatments (including a "no treatment" option) prior to performing any procedure.
- Use of proper patient selection criteria is critical to clinical success:
 - Patients shall be adults, from the age of 18 years old.
 - Ethnicity doesn't have any impact on the use of these medical devices.
 - Due to the potential effect of the anaesthesia, do not use on pregnant women.
 - The patient's current health status, medical and dental history have to be evaluated taking into account the contraindications, warnings and precautionary information.
 - Determine the location of all the anatomical features to be avoided prior to initiating any implant procedure.
 - Patients with hyperfunction (bruxism) may be at increased risk for implant failure or fracture.
 - Adequacy of bone should be determined by X-rays, palpation, probing, and thorough visual inspection of the proposed implant site.
 - Poor bone quality, poor patient oral hygiene habits, tobacco use, certain concomitant medications, and systemic disease (diabetes, etc.) may contribute to lack of osseointegration and subsequent implant failure.
- DURA-VIT Dental Implants are sterile and single-use medical devices. Do not sterilise and re-use. Indeed, any re-use and re-sterilization of the medical devices are strictly prohibited as the functional performance and safety of the device expected by the manufacturer cannot be guaranteed by exposing the patient to the risk of loss of device performance and/or serious risks of contamination and cross-infections.

Proper use of surgical instrumentation is critical to clinical success:

- The surgical instrumentation has been specifically designed for the dental implant lines by B. & B. Dental.
- Any use and re-use of surgical instruments without an accurate and in-depth reading of their specific leaflets does not guarantee the functional performance and safety of the devices provided by the manufacturer, thus exposing the patient to the risks of non-mechanical performance of the product and / or serious risks of contamination.
- All surgical instruments must be checked before each use regarding their state of cleanliness and sterilisation to ensure that they are always in operating condition. Do not use devices that show anomalies as they could compromise the safety and health of the patient.
- Inspect all surgical instruments before use. Using dull or worn instruments may cause implants to fracture. Use of excessively worn drivers can cause lack of implant advancement or their locking onto the implant head.
- Ensure there is a solid connection between any driver, handpiece or torque ratchet and implant.
- Take adequate precautions to avoid that the patient inhales and/or ingests small components.
- Care must be taken when using instrumentation near other teeth to avoid damaging tooth structure.
- During implant placement and fastening of the screws do not exceed the torques indicated in the paragraph "**CAUTION: Screwing instructions**".
- Excessive torque in hard bone may result in implant fracture or over-compression of bone, which could lead to necrosis.
- Patients with nocturnal bruxism who have removable prostheses should sleep with the prosthesis in place or wear a protective mouth-guard covering the abutment head(s).
- Electrosurgery is not recommended around metallic implants or the abutments, due to risk of electrical and/or heat conductivity.
- The pure titanium and titanium alloy used to manufacture DURA-VIT implants are not magnetic and thus should be acceptable for MRI procedures. However, it is important that patients advise medical professionals about the presence of dental implants prior to imaging to ensure compatibility with the MRI equipment and to address potential artifacts, depending on the area of imaging.
- Medical professionals, who performed the surgery, must apply the traceability label, present in the cardboard box, in the dedicated field on the implant card and notify the patient.
- Instruct the patient to follow a routine post-surgery regimen as reported on the implant card, according to the paragraph "**Post-surgery course**".
- Loading a dental implant beyond its functional capacity can result in a fractured implant, loss of bone, failure to osseointegrate, and/or subsequent loss of osseointegration.
- Do not probe excessively around the implant, as this can lead to soft tissue damage and peri-implantitis.
- Complications

- Unsuccessful osseointegration is evidenced by infection, implant mobility or bone loss.
- Remove any failed implant, and any surrounding visible granulation tissue, from the implant site as soon as possible.

Storage and use

- All devices must be stored carefully in a clean and dry environment to keep the package intact.
- Devices should be stored at room temperature.
- Refer to individual product labels for special storage and handling conditions.
- Sterile products should be used before the expiration date printed on the label.
- Devices must not be used whenever the packaging is not intact.
- Contents are considered sterile unless the package is opened or damaged. Do not use products labelled "STERILE" if packaging has been damaged or opened prior to intended use.

Disposal

Dispose of in accordance with the laws in force as special hospital waste. In particular, if the devices are contaminated with blood or fluids, it is necessary to use the appropriate containers and gloves to avoid any direct contact.

Surgical Protocol for Implant Placement

Patient selection

Careful case planning and appropriate patient selection criteria are critical and extremely important for successful clinical procedures. This should include assessment and review of the patient's current health status, medical and dental history, and the state and ability to maintain proper oral hygiene. Anatomical acceptability must be evaluated, implying that the adequacy of bone shall be determined by X-rays, palpations, probing, and thorough visual inspection of the proposed implant site. It is essential also to determine the location of all the anatomical structures that must be avoided before beginning any type of implant procedure. In particular, clinical and radiographic examination of ridge thickness and soft tissue condition is required. It is always recommended to perform a cone-beam CT scan to have precise anatomical references of thickness, height, density of bone and tissues.

Implant selection

The choice of the size of the implant is important for long-term success. It is better to use the maximum diameter and length possible, for greater stability of the superimposed prosthesis. Use X-rays and radiographic sizing transparencies to determine the best implant type, diameter and length. Evaluate the bone density and soft tissue depth to determine the type of implant to use.

The measurements can be made directly on the panoramic X-ray, superimposing the transparencies on it, enlarged by 25% compared to the real image. The edges of the crest must be adequately palpated in order to estimate an insertion angle that must reach a good parallelism with the other implants.

In case of cone-beam CT scan, all the necessary anatomical references are immediately available for making the evaluation.

	Bone Density				Soft Tissue Depth		Position			Implant placement	
	D1	D2	D3	D4	< 2 mm	≥ 2 mm	Frontal	Lateral	Molar	Under crest	Crest
DURA-VIT Implant line											
SLIM Ø 3,0	✓	✓	✓	✓	✓	✓	✓	NR	NR	✓	✓
SLIM Ø 3,4	✓	✓	NR	NR	✓	✓	✓	NR	NR	✓	✓
3P	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
EVOLUTION	NR	NR	✓	✓	✓	✓	✓	✓	✓	✓	✓
WIDE	NR	✓	✓	✓	✓	✓	NR	✓	✓	✓	✓
PTERYGO	✓	NR	✓	✓	✓	✓	NR	NR	✓	✓	✓

*D1 = Very Dense Bone

*D4 = Very Soft Bone

NR = Not Recommended

The appropriate implant diameter should be completely surrounded by at least 1.0 mm of bone and the length should engage bone for the entire threaded portion of the implant. A 1.0 mm subcrestal implant placement is also recommended.

Patient preparation

Prior to preparing the patient for surgery, take an occlusal bite registration. Prepare the patient, following standard aseptic guidelines. Administer sufficient anaesthetic by local infiltration. In the mandible, inferior alveolar nerve block anaesthesia is not recommended in order to allow the patient to retain sensation and give feedback when necessary.

Surgical procedure

As with any other surgical action, it is important for the implant procedure to be performed in an environment which is as sterile as possible. The surgical instrumentation is supplied non-sterile. It shall be disinfected, cleaned and sterilized before each use according to the specific IFU.

General instructions for preparing the implant site

An efficient and atraumatic implant site preparation is created through a procedure relying on a gradual drilling technique. The whole stage of bone tissue drilling must be performed under an abundant external irrigation with saline solution or, preferably, sterile double-distilled water. Furthermore, drilling must be intermittent both to avoid bone to heat up and to create a pumping effect that will help effective removal of bone tissue.

STANDARD SURGICAL PROTOCOL (lines DURA-VIT SLIM, 3P, EV, WIDE and PTERYGO)

Preparation of the site in open flap surgery

Make a full thickness mesiodistal incision along the surface of the crest. Using a periosteal elevator, leave the surface exposed to the alveolar bone only as necessary in order to have adequate vision and adequate working space for the surgical procedure. Pointed crest or other bone irregularities must be removed using the proper surgical instrument, such as a leveller. A minimum distance of 2 and 6 mm (from cut to cut) must be maintained between two implants and/or implant and the adjacent natural teeth.

The implant site must be prepared in a visible, clear field, so that the operator can prepare it with precision in every phase, thus allowing for proper positioning of the implant.

Site preparation in bone D1 and D2

1. FIRST DRILL (guide drill or lance drill) allows for the direct perforation of the bone cortex at the centre of the crest. Operating speed: 800-1000 RPM.
2. SITE LENGTH PREPARATION: The depth drill, size 2.1 mm, equipped with marked stops indicating the depth, prepares the implant site of the selected length. Operating speed: 800-1000 RPM.
3. SITE WIDTH PREPARATION: The shaping drills for the diameter of the selected implant are calibrated with the core of the implant. They have to be used sequentially by increasing diameter, from 3.0 mm up to the diameter of the implant to be inserted, in order to obtain the maximum precision with the minimum possible heating. Removable stops are available for every depth and widening drill to facilitate the implant site preparation. Operating speed 100-750 RPM.
4. PREPARATION OF THE SHOULDER in cortical bone: the preparation phase ends with the use of countersink drills of the correspondent diameters.

Site preparation in bone D3 and D4

1. FIRST DRILL (guide drill or lance drill) allows for the direct perforation of the bone cortex at the centre of the crest. Operating speed: 600-800 RPM.
2. SITE LENGTH PREPARATION: The depth drill, size 2.1 mm, equipped with marked stops indicating the depth, prepares the implant site of the selected length. Operating speed: 600-800 RPM.
3. SITE WIDTH PREPARATION: The use of compactor-expanders is required. The shaping compactors for the diameter of the selected implant are calibrated with the core of the implant. They have to be used sequentially by increasing diameter, from 3.0 mm up to the diameter of the implant to be inserted, in order to obtain the maximum precision with the minimum possible heating. For facilitating the implant site preparation, the same stops available for drills have to be used for every compactor. Operating torque: maximum 70 Ncm.

GUIDED SURGERY PROTOCOL FOR THE DURA-VIT IMPLANT LINES SLIM, 3P, EV AND PTERYGO

All the surgical instruments for the guided surgery procedure have a guide of 9.0 mm of height and two different diameters, 4.2 and 5.5 mm, which allow the guided implant site preparation on the base of the project previously developed on the CT scan processed according to the specific protocol by a dedicated software.

Preparation of the site

The guided surgery protocol allows the flapless procedure for the implant site preparation.

Therefore, when adherent mucosa is abundant on the proposed implant site, it is recommended to use the guided tissue punch to remove tissue and gain access to the bone.

Otherwise, if it is thin, the gingiva and bone can be drilled directly with the guided lance drill.

In guided surgery it is necessary to use a surgical guide, in which the proper sleeves must be fixed, according to the diameter of the implant to be placed. Indeed, they are available in two different diameters for the same length (5.0 mm): 4.2 mm and 5.5 mm. 4.2 mm sleeves are suitable for implants with a diameter from 3.0 mm to 4.0 mm. The wider ones are for implants with a diameter of 4.5 or 5.0 mm and they require the use of the converter attached to the drills for the drill stability and to avoid any movement.

Site preparation in bone D1 and D2

1. FIRST DRILL (guide drill or lance drill) allows for the direct perforation of the bone cortex at the centre of the crest. Operating speed: 800-1000 RPM.
2. PREPARATION OF THE SHOULDER. The preparation phase starts with the use of countersink drills, from size 3.0 mm up to the implant diameter. Operating speed: 800-1000 RPM.
3. SITE WIDTH PREPARATION. The shaping drills for the diameter of the selected implant are calibrated with the core of the screw. They have to be used sequentially by increasing the diameter, from 2.1 mm up to the diameter of the implant to be inserted, in order to obtain the maximum precision with the minimum possible heating. Operating speed: 800-1000 RPM.
4. SITE LENGTH PREPARATION. Use the drill of the same length of the selected implant and bring it to the end of its stroke, making sure to use cooling systems to avoid excessive overheating. The drills prepare an osteotomy increased by 0.5 mm in relation to the length.

Site preparation in bone D3 and D4

1. FIRST DRILL (guide drill or lance drill) allows for the direct perforation of the bone cortex at the centre of the crest. Operating speed: 600-800 RPM.
2. PREPARATION OF THE SHOULDER. The preparation phase starts with the use of countersink drills, from size 3.0 up to the implant diameter minus 0.5 mm, since in case of spongy bone the implant site has to be used underprepared by 0.5 mm. Operating speeds: Operating speed: 350-600 RPM.
3. SITE WIDTH PREPARATION. The shaping drills for the diameter of the selected implant are calibrated with the core of the implant. They have to be used sequentially by increasing the diameter, from 2.1 mm up to the diameter of the implant to be inserted minus 0.5 mm, in order to obtain an underprepared implant site with the maximum precision and the minimum possible heating. Please note that the osteotomy corresponds to the implant diameter, if the compactor-expanders are employed. Operating speeds: 350-600 RPM.
4. SITE LENGTH PREPARATION. Use the drill or the compactor-expander of the same length of the selected implant and bring it to the end of its stroke, making sure to use cooling systems to avoid excessive overheating. The drills prepare an osteotomy increased by 0.5 mm in relation to the length.

Treatment and sterilisation

The implants are supplied in sterile packages, and they must be placed with cleaned and sterile surgical instruments. Given the nature of surfaces treatment (sandblasting or double etching), DURA-VIT Dental Implant lines have to be handled with titanium instruments and/or sterile, powder-free gloves to avoid contaminating the product or interfering with its biocompatibility.

Implant placement

1. Irrigate the entire implant site with saline solution or, preferably, sterile double-distilled water.
2. Secure the implant with the implant driver for contra-angle or for torque ratchet, extract it and then bring it into the prepared site. Rotate the implant clockwise until it is fixed.
3. The cover screw and transmucosal screw, if provided, are included in the cap of the implant, and both must be tightly screwed onto the implant before suturing the eventual flap.
4. An X-ray check of the position of the implant is recommended before the potential closure of the flap.
5. The eventual mucoperiosteal flap must be delicately repositioned for optimal adaptation of the tissue, and then sutured.

WARNING: Screwing instructions

Please respect the insertion forces for the placement of the implant and of the screws indicated here below:

- Maximum 35 Ncm with implant drivers for contra-angle.
- Maximum 70 Ncm with implant drivers for torque ratchet.
- Only for SLIM implants it is necessary not to exceed a force of 45 Ncm with implant drivers for torque ratchet.
- Maximum 10 Ncm with cover screw.
- Maximum 20 Ncm with transmucosal and healing screws.

Post-surgery course

The patient must be taught to follow the postoperative course. The patient should be informed on how to behave immediately after surgery and in the subsequent days, which includes cold compresses, no physical exertion and/or sauna for the first 24 hours, no alcohol, no nicotine, no tea or coffee in the first days. The patient must be cautioned not to chew in the implant area. Moreover, s/he should be advised that poor oral hygiene, tobacco smoking and systemic and general diseases (diabetes, rheumatoid arthritis, etc.) can contribute to poor bone integration and a consequent failure of the implant.

The surgeon could evaluate the possible administration of medicines pre- and post-surgery.

The stitches must be removed after one week.

It is recommended that every removable prosthesis which remains near the implant site be adequately raised and aligned using a soft material for regeneration.

Later, the patient shall maintain an appropriate and continuous oral hygiene, have regular check-ups and professional cleaning, and contact their dentist if s/he experiences tender points, note loss of retention or bone resorption.

Patients shall advise medical professionals about the presence of the implanted medical devices prior to MRI. In most cases, these medical devices will not interfere with MRI.

Customer information

No person is authorised to provide any information which deviates from the information provided in this instruction sheet.

For further information and technical documents, please consult our website www.bebdental.it, especially the "download" page.

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED, through the following website: <https://ec.europa.eu/tools/eudamed>. The SSCP is linked to the following BASIC UDI-DI:

- **805228249BIPHASICIMPLAN6P (biphasic implants);**
- **805228249MONOPHASICIMPLANB6 (monophasic implants);**
- **805228249MINIIMPLANTS DP (mini implants);**
- **805228249TEMPORARYABTCL (titanium temporary abutments);**
- **805228249PTEMPORARYABT7G (castable and peek temporary);**
- **805228249FINALABUTMENTXU (final abutments);**
- **805228249SPHERICALABTBZ (spherical abutments);**
- **805228249HEALINGDEVICESWA (healing screws, cover and transmucosal screws).**

In the event that serious or suspected accidents occur following the use of our medical devices, you shall report them immediately to the Company and to your National Competent Authority.

Table 1

CE Mark with the number of the Notified Body TÜV Rheinland	
CE Mark of a medical device	
Date of Manufacture	

Use by	
Lot Code	
Device Number	
Manufacturer	
Consult Instruction of Use	 https://www.bebdental.it/pro/en/instructions/
Do Not Re-use	
Caution	
Sterilised using irradiation Single sterile barrier system with protective packaging outside	
Do not use if package is damaged or open	
Unique Device Identifier	
Do not resterilize	
Non-Sterile Device	
Medical device	
Keep away from the sunlight	

Keep dry	
----------	---

Table 2

Patient name or patient ID	
Date of implantation	
Name and Address of the implanting healthcare institution/provider	
Information website for patient	
Medical device	

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