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
INSTRUCTIONS FOR USE

General Information & Intended Use

Mini dental implants of the line DURA-VIT MINI are implantable medical devices used in case of total or partial mandibular and/or maxillary edentulism and/or dental agenesis. They get in contact with bone and mucous tissue to replace one or more missing natural teeth, up to the reconstruction of the entire dental arch. They are a support for a fixed prosthesis or an anchoring of removable prosthesis. The expected benefit is both functional, to rehabilitate correct and complete chewing, as well as to restore the ability to speak, and aesthetic, such as smiling.

Mini dental implants of the line DURA-VIT MINI utilize a self-tapping, threaded screw design and employs minimally invasive surgical intervention. The system includes small diameter dental implants with components enabling fixed and removable restorations, as well as metal housings, abutment copings, surgical and prosthetic instrumentation, and laboratory components.

Mini dental implants of the line DURA-VIT MINI are available with the following specifications:

- Diameters: 2.0 mm, 2.4 mm, 2.5 mm
- Intra-osseous lengths: 10 mm, 13 mm, 15 mm
- Abutments: O-ball or square head
- Available with 2.5 mm transgingival collar to accommodate various soft tissue depths.

Indications

- Long-term partial and full denture stabilization
- Long-term fixation of bridges

Material in use

Mini dental implants are realized in titanium alloy Ti6Al4V (grade 5).

Materials results compliant with harmonized standards.

The allergy to titanium is a very rare but possible event, so it is always necessary to check in advance with patients, who do not have this type of allergy.

Warning

- Under law, only surgeons and certified dentists are allowed to use implants in surgery. Any use by other persons is strictly prohibited. Therefore, DURA-VIT MINI Implant System must be handled and placed only by qualified and specially trained personnel.
- B. & B. Dental assumes no responsibility for any possible damages to the patient or any serious health hazards caused through improper use, or any use by unqualified or inadequately trained personnel, that is not in accordance with this instruction manual.
- The manufacturer excludes any responsibility for any damage to the patient from the use of medical devices not manufactured and not supplied by B. & B. Dental itself.
- Medical devices of this IFU must be placed in a hospital or clinic environment.
- Please check "Table 1" for a better understanding of the symbols on medical devices labels.
- Please check "Table 2" for a better understanding of the symbols on the implant card.
- Potential adverse events associated with the use of this product may include, but are not limited to, the following.
 - Fractured and broken implant
 - Failure of or inadequate osseointegration
 - Unintended perforation of bony and/or soft tissue, such as inferior alveolar canal, maxillary sinus, and lingual artery.
 - Fracture of labial or lingual plate
 - Partial nerve injury or irreversible damage (paraesthesia/hyperesthesia/anaesthesia)
 - Loss of bone due to concomitant events, local or systemic infection, peri-implantitis, loss of soft tissue volume

Proper selection of patient candidates, implant placement and post-operative care, combined with correct product use may help to mitigate these risks.

General contraindications

Dental implants shall not be placed in case of noncompliant patient and if there is any intraoral contraindication, such as anatomical impairments, poor quality bone tissue or insufficient alveolar bone quality and quantity to achieve initial stability, or if there is clinical or radiographic evidence of pathology.

Absolute contraindications

Known allergies to titanium, recent myocardial infarction, valvular prosthesis, severe renal disorder, severe liver disorders, osteomalacia, treatment resistant diabetes, recent high-dosage radiation therapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.

Relative contraindications

Chemotherapy, moderate renal disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged cortico-steroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Precautionary Information

- Since the surgical techniques are highly specialized and complex procedures, it is strongly recommended that practitioners shall complete a training course to learn established techniques of placement in oral implantology. Failure to seek proper training or follow Instructions for Use (IFU) may impact clinical success of implant, such as failure of the implant, loss of bone, and postoperative complications, that could lead to pain for the patient and put him under serious risks his health.
- Surgical placement and restoration of dental implants requires intricate procedures that may involve risks. It is the responsibility of the clinician to inform the patient of any risk and alternative treatments (including a "no treatment" option) prior to performing any procedure.
- It is recommended that clinicians who place mini implants be capable of performing flapped procedures.
- Use of proper patient selection criteria are critical to clinical success:
 - Patients shall be adults, from the age of 18 years old.
 - Ethnicity doesn't have any impact on the use of these medical devices.
 - Due to the potential effect of the anaesthesia, do not use on pregnant women.
 - The patient's current health and psychological status, medical and dental history have to be evaluated taking into account the contraindications, warnings and precautionary information.
 - Determine the location of all avoidable anatomical features prior to initiating any implant procedure.
 - Patients with hyperfunction (bruxism) may be at increased risk for implant failure or fracture.
 - Adequacy of bone should be determined by radiographs, palpation, probing, and thorough visual inspection of the proposed implant site.
 - Insufficient availability of bone, poor bone quality, poor patient oral hygiene habits, tobacco use, certain concomitant medications, and systemic disease (diabetes, etc.) may contribute to lack of osseointegration and subsequent implant failure.
- Mini implants DURA-VIT MINI are sterile and single-use medical devices. Do not sterilize and re-use. Doing so may result in product failure and potential patient injury. Indeed, any re-use and re-sterilization of the medical devices are strictly prohibited as the functional performance and safety of the device expected by the manufacturer cannot be guaranteed by exposing the patient to the risk of loss of device performance and / or serious risks of contamination and cross-infections.

Proper use of surgical instrumentation is critical to clinical success:

- The surgical instrumentation has been specifically designed for the dental implant line by B. & B. Dental.
- Any use and re-use of surgical instruments without an accurate and in-depth reading of its specific leaflet does not guarantee the functional performance and safety of the devices provided by the manufacturer, exposing the patient to the risks of non-mechanical performance of the product and / or serious risks of contamination.
- All surgical instruments must be checked before each use regarding their state of cleanliness and sterilization to ensure that they are always in an operating condition. Do not use devices that show anomalies as they could compromise the safety and health of the patient.
- Inspect all surgical instruments prior to use. Using dull or worn instruments may cause implants to fracture. Use of excessively worn drivers can cause lack of implant advancement or their locking onto the implant head, or stripping of the abutment.
- Ensure there is a solid connection between any driver, handpiece or torque ratchet and implant.
- Take adequate precautions to avoid patient aspiration and/or ingestion of small components.
- Care must be taken when using instrumentation near other teeth to avoid damaging tooth structure.
- Do not exceed 45 Ncm of torque during implant placement. Excessive torque may result in implant fracture or over compression of bone, which could lead to necrosis.
- Excessive torque in hard bone may result in implant fracture or over compression of bone, which could lead to necrosis.

- Patients with nocturnal bruxism who have removable prostheses should sleep with the prosthesis in place or wear a protective mouth-guard covering the abutment head(s).
- Electro surgery is not recommended around metallic implants or the abutments, due to risk of electrical and/or heat conductivity.
- The titanium alloy used to manufacture mini-implant DURA-VIT MINI are not magnetic and thus should be acceptable for MRI procedures. However, it is important that patients advise medical professionals about the presence of dental implants prior to imaging to ensure compatibility with the MRI equipment and to address potential artifacts, depending on the area of imaging.
- Medical professionals, who performed the surgery, must apply the traceability label, present in the cardboard box, in the dedicated field on the implant card and notify the patient.
- Instruct the patient to follow a routine post-surgery regimen as reported on the implant card, according to the paragraphs "**Post-surgery course**" and "**Patient Home Care & Maintenance Instructions**".
- Loading a dental implant beyond its functional capacity can result in a fractured implant, loss of bone, failure to osseointegrate, and/or subsequent loss of osseointegration.
- Do not probe excessively around implant, as this can lead to soft tissue damage and peri-implantitis.
- Complications
 - Unsuccessful osseointegration is evidenced by infection, implant mobility or bone loss.
 - Remove any failed implant, and any visible granulation tissue, from the implant site as soon as possible.

Storage and use

- All devices must be stored carefully in a clean and dry environment to keep the package intact.
- Devices should be stored at room temperature.
- Refer to individual product labels for special storage and handling conditions.
- All titanium implants and parts must be handled via titanium locking pliers and/ or sterile, powder-free gloves to avoid contaminating the product or interfering with its biocompatibility.
- Sterile products should be used before the expiration date printed on the label.
- Devices must not be used whenever the packaging is not intact.
- Contents are considered sterile unless the package is opened or damaged. Do not use products labelled "STERILE", if packaging has been damaged or opened prior to intended use.
- Reusable instruments are provided non-sterile and they shall be cleaned, disinfected and sterilized according to the IFU.

Disposal

Dispose of in accordance with the laws in force as special hospital waste. In particular, if the devices are contaminated with blood or fluids, it is necessary to use the appropriate containers and gloves to avoid any direct contact.

Surgical Protocol for Implant Placement

Patient Selection:

Careful case planning and appropriate patient selection criteria are necessary for successful procedures of DURA-VIT MINI. This should include a review of the patient’s current health status, medical and dental history, and the ability to maintain proper oral hygiene. In particular, physical and radiographic examination of ridge thickness and soft tissue condition is required. Ridge mapping calipers under local anaesthetic may be used to create a map of bone thickness along the arch and location of key anatomical structures. While most cases can be done without cutting a soft tissue flap, clinical judgment should be used to determine if the case warrants a flap procedure. It is recommended that clinicians who place mini-implants be capable of performing flapped procedures.

Implant Selection:

Use radiographs and radiograph sizing transparencies to determine the best implant type, diameter, and length. Evaluate the bone density, soft tissue depth, and buccolingual width to determine the type of implant to use.

Implant type	Bone Density				Soft Tissue Depth		Buccolingual Width		
	D1	D2	D3	D4	< 2 mm	≥ 2 mm	≤ 4 mm	>4 mm	≥ 5 mm
2.0 mm with Collar	✓	✓	NR	NR	NR	✓	NR	✓	✓
2.4 mm with Collar	NR	✓	✓	NR	NR	✓	NR	NR	✓

2.5 mm with Collar	NR	✓	✓	NR	NR	✓	NR	NR	✓
--------------------	----	---	---	----	----	---	----	----	---

*D1 = Very Dense Bone
 *D4 = Very Soft Bone
 NR = Not Recommended

The appropriate implant diameter should be completely surrounded by at least 1.0 mm of bone and the length should engage bone for the entire threaded portion of the implant.

Determination of Implant Sites in the Mandible:

The implants should be placed at least 7 millimetres anterior to the mental foramina. A minimum of four mini implants should be placed at least 5 millimetres apart to accommodate the metal housings. Mark the proposed site for each implant with a bleeding point or intraoral marker. Implant fracture can occur due to excessive load. Treatment planning should consider the need for additional implants based on occlusal forces and occlusal load.

Determination of Implant Sites in the Maxilla:

If placing in the maxilla, evaluate position of sinuses in relation to the ridge. A minimum of six mini implants are recommended in the maxilla, placed at least 5 millimetres apart to allow space for the metal housings. Mark the proposed site for each implant with a bleeding point or intraoral marker. Implant fracture can occur due to excessive load. Treatment planning should consider the need for additional implants based on occlusal forces and occlusal load.

Patient Preparation:

Prior to prepping the patient for surgery, take an occlusal bite registration. Prep the patient, following standard aseptic guidelines. Administer sufficient infiltration anaesthetic. In the mandible, inferior alveolar nerve block anaesthesia is not recommended in order to allow the patient to retain sensation and give feedback in instances of implant encroachment. In a flapless procedure when adherent mucosa is present at the proposed implant site, it is recommended to use a 1.5 millimetre tissue punch to remove tissue and gain access to the bone.

Implant Placement:

Always place the two most distally-sited implants first, and work toward the midline.

Use the 1.1 mm surgical drill, with adequate external, sterile irrigation, at the first implant site to create a pilot hole by pumping the drill slightly up and down at a speed of 1200 to 1500 rpm until the cortical plate is penetrated. Continue to carefully advance the drill to a depth of one third to one half the length of the threaded portion of the selected implant. Do not drill beyond one half depth in D2-D3 bone because the threaded implant is self-tapping and must use the non-drilled bone to provide immediate fixation.

NOTE: In some dense bone (D1) cases, additional drilling depth and/or use of one or more drills of slightly larger diameter may be required to seat the implant without exceeding 45 Ncm of torque. See Troubleshooting instructions below. Grasping the attached blue vial cap finger driver, remove the implant from the vial taking care not to touch the sterilized implant body. You may use the cap to carry the implant to the site. Carefully align the implant to the pilot hole and insert the tip. Turn the implant clockwise with even downward pressure. The implant threads will engage bone. Align the implant as vertically as possible. It is possible to slightly modify the alignment of the implant to the ridge, thus improving parallelism with subsequently placed implants. Throughout the placement procedure, the implant should be advanced slowly into the bone to allow time for bone relaxation and to minimize thermal effects. Using the blue cap, make the first turns until the cap comes off independently. Then replace it with the wing screwdriver and turn the implant deeper into the bone. Be careful not to exceed the torque of 45 Ncm at the implant. It is possible to bring the implant into its end position with the wing screwdriver. However, the graduated torque ratchet should be used for the final rotations to have an accurate record of implant stability at final placement and to ensure that the recommended maximum torque of 45 Ncm is not exceeded.

NOTE: If the implant turns easily, there may not be sufficient bone density to securely place or immediately load the implant. See Troubleshooting Instructions below. Replace the Winged Thumb Wrench with the graduated torque wrench assembly (see Torque Wrench IFU for assembly instructions). Slowly and carefully make small turns to advance the implant to its final position with the abutment head protruding from the gingival tissue and all surface-roughened threads embedded in bone. **At final placement, a minimum of 35 Ncm torque must be achieved to allow immediate loading of the implant. 45 Ncm torque should not be exceeded as this could cause bone damage.**

Once the implant is in its final position, proceed to the next implant site and repeat the above process. Align subsequently placed implants as parallel as possible. The O-ring denture fixation system allows up to a 30 degree range of misalignment. However, the more parallel the implants are to each other, the better the connection and fit, resulting in less wear on the system.

Troubleshooting

In cases where the implant turns easily or torsional forces approaching 45 Ncm would be required to fully seat the implant, consider the alternatives below:

Fully seated implant sing torque <35 Ncm	Incompletely Seated Implant with Torque Approaching 45 Ncm
<ul style="list-style-type: none"> ▪ Use Soft Reline protocol rather than immediate loading protocol. ▪ Replace with larger diameter and/ or longer implant. ▪ Move implant to alternative site. 	<ul style="list-style-type: none"> ▪ Extend drill depth of pilot hole (not to exceed length of implant). ▪ Use a drill larger than 1.3 mm of diameter to enlarge the osteotomy. ▪ Replace with smaller diameter and/or shorter implant. ▪ Move implant to alternative site.

NOTE: Use caution when removing an implant to avoid implant fracture.

NOTE: Implant placement into extremely soft (D4) bone is not recommended.

Metal Housing Selection

Select the appropriate housing per the table below.

NOTE: Ensure a minimum of 1 mm of acrylic surrounds the housings in the denture.

	MD-3004/1	MD-3004	MD-3004/2
Height	3.5 mm	3.3 mm	3.1 mm
Diameter	4.7 mm	4.3 mm	4.0 mm
Angulation Tolerance	Up to 30°	Up to 15°	0 - 5°

To replace worn O-rings, purchase O-rings separately. See 'O-ring Replacement' instructions below.

Chairside Prosthetic Protocol for Housing Attachment

In mandibular cases, once all implants have been placed and exhibit a minimum of 35 Ncm of stability, the denture can be immediately adapted to the implant abutments in a hard pick-up procedure. In the maxilla, use of secure soft reline is recommended for the first 6 months. See Chairside Prosthetic Protocol for Soft Reline Instructions below.

NOTE: If a flapped procedure was utilized, protect sutures using a rubber dam prior to performing the hard pick-up procedure. First, place a housing on each implant. Determine their locations in the fitting surface of the denture, using pressure indicating paste, bite registration material, or a similar method. Then create individual holes or a trough to accommodate the housings with a passive fit to the denture. You may also choose to relieve the denture border a few millimeters for added patient comfort. For each implant, cut a length of blockout shim to provide a seal between the tissue and the bottom of the housing. This is to block out any undercuts. Place the shim over the implant and replace the housing, then reconfirm the passive fit.

NOTE: Ensure that the shim is not preventing the housing from fitting all the way onto the implant. Remove the denture. Clean and dry the modified denture surface. Petroleum jelly can be applied to desired areas to avoid unintended bonding of the reline material. Next, follow the IFU for the secure hard pick-up material (or other chairside reline material), and apply the material into the relieved space of the denture and directly onto the top of the housings intraorally. Use the bite registration to properly align the denture to the preoperative occlusal relationship. Carefully insert the denture over the housings and have the patient apply normal biting pressure for 7 to 9 minutes until the secure hard pick-up has firmly set. (Check the manufacturer IFUs for cure times of other reline materials.) This creates a firm fixation of the housings to the denture at the correct position for each abutment. Remove the denture with the now-embedded housings. Remove the blockout shims. Trim, fill if necessary, and polish the tissue-contact surfaces of the denture to achieve a comfortable, firm fit. For maxillary dentures, the palate can be removed, if desired, after osseointegration is achieved. After final fitting, instruct the patient to leave the denture in place for 48 hours to prevent tissue overgrowth on the implant abutments. Patient should be advised to eat soft foods for the first few days.

Chairside Prosthetic Protocol for Soft Reline

In the mandible, if the torque at final positioning of one or more implants is below 35 Ncm, a soft reline method for retaining the denture is recommended for 4-6 months or until osseointegration has occurred. **In the maxilla, use of a soft reline is recommended for the first 6 months in all cases.** See Chairside Prosthetic Protocol for Soft Reline Instructions below.

NOTE: If a flapped procedure was utilized, protect sutures using a rubber dam prior to performing a soft reline procedure. Relieve the denture base at least 1 mm and relieve the tissue contact surface to accommodate the

prosthetic heads of each implant. Roughen the surface with an acrylic bur and degrease the surface with isopropyl alcohol. Next, follow the IFU for the secure soft reline (or other chairside soft reline material), and apply the material into the relieved space of the denture. Use the bite registration to properly align the denture to the preoperative occlusal relationship. Carefully insert the denture over the implants and have the patient apply normal biting pressure for 10 minutes until the secure soft reline has firmly set. (Check the manufacturer IFUs for cure times of other soft reline materials.) This will create abutment orifices aligned with the implants in the soft reline material. Remove the denture and trim excess material. Follow finishing instructions provided by the soft reline material manufacturer.

For maxillary dentures, do not remove the denture during this stage

Instruct the patient to keep the denture in place for the first 48 hours after placement to prevent soft tissue overgrowth. It is recommended to wait several weeks before eating hard foods to allow bone in-growth around the implant. After four to six months when the implants have osseointegrated, the soft liner can be replaced with metal housings following the Chairside Prosthetic Protocol for Housing Attachment.

Indirect prosthetic protocol for housing attachment or single unit crown

Once all implants have been placed and exhibit a minimum of 35 Ncm of stability, snap the o-ball impression copings directly onto the prosthetic heads of the mini implants. Copings can be modified as necessary but do not remove external retentive features.

NOTE: Soft tissue may prevent full engagement of the coping on implants seated too deeply into soft tissue. In such cases, it is recommended to take an impression of the prosthetic head of the implant without impression copings applied. Use an impression material with good flow properties and adequate rigidity. Standard crown and bridge impression techniques are used to pick up the impression copings, recording each implant's position. Refer to the impression material IFU.

NOTE: If a flapped procedure was utilized, protect sutures using a rubber dam prior to taking the impression. Once the impression has fully set, carefully remove the tray from the patient's mouth and confirm all impression copings have been captured accurately in the impression. Select the appropriate lab analogs - refer to product catalog. Insert the analogs into the copings. This step can be accomplished in the clinic or at the dental laboratory. Use standard stone model fabrication techniques to form the model.

Post-surgery course

The patient must be taught to follow the postoperative course. The patient should be informed on how to behave immediately after surgery and in the subsequent days, which includes cold compresses, no physical exertion and/or sauna for the first 24 hours, no alcohol, no nicotine, no tea or coffee in the first days). The patient must be cautioned not to chew in the area where the membrane was inserted. Moreover, he should be advised that poor oral hygiene, tobacco smoking and systemic and general diseases (diabetes, rheumatoid arthritis, etc.) can contribute to poor bone integration and a consequent failure of the implant.

The surgeon could evaluate the possible administration of medicines pre- and post-surgery.

The stitches must be removed after one week.

It is recommended that every removable prosthesis which remains near the implant site be adequately raised and aligned using a soft material per regeneration.

Patient Home Care & Maintenance Instructions

Inform the patient that they need to treat their new implants very much like teeth, which includes the following.

1. Regularly brush the implants to keep them free of plaque and debris.
2. Maintain an appropriate and continuous oral hygiene.
3. Clean and sanitize dentures regularly. We recommend brushing and soaking overnight using a solution of equal parts vinegar and water. Avoid common cleaning chemicals including bleach, hydrogen peroxide, and solvents of any kind.
4. Return for regular check-ups and professional cleaning.
5. Contact your dentist if you experience sore spots, notice loss of retention, dentures rock back and forth, or an O-ring is lost or damaged.
6. Advise medical professionals about the presence dental implants prior to MRI imaging. In most cases, titanium dental implants will not interfere with MRIs.

Maintenance Instructions for the surgeon

- During regular check-ups and cleaning the fit and function of the denture shall be verified, replace worn O-rings and reline as necessary.
- In cases of bone resorption, adjustment of the denture is required to ensure implants are not subjected to overloading. The DURA-VIT MINI Implant System is intended to be a tissue supported, implant retained system not an implant supported system.

O-ring Replacement

The connection between the O-Ball abutment on the implants and the metal housings may lose retention over time due to O-ring wear.

1. Remove O-ring from the housings using a dental explorer or similar instrument.
2. Insert a new O-ring by aligning the O-ring over the opening of the housing and pushing into place with a ball burnisher or another rounded instrument.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED, through the following website: <https://ec.europa.eu/tools/eudamed>. The SSCP is linked to the following BASIC UDI-DI:

- 805228249BIPHASICIMPLAN6P (biphasic implants);
- 805228249MONOPHASICIMPLANB6 (monophasic implants);
- 805228249MINIIMPLANTS DP (mini implants);
- 805228249TEMPORARYABTCL (titanium temporary abutments);
- 805228249PTEMPORARYABT7G (castable and peek temporary);
- 805228249FINALABUTMENTXU (final abutments);
- 805228249SPHERICALABTBZ (spherical abutments);
- 805228249HEALINGDEVICESWA (healing screws, cover and transmucosal screws).

In the event that serious or suspected accidents occur following the use of our medical devices, you shall report them immediately to the Company and to your National Competent Authority.

Table 1

CE Mark with the number of the Notified Body TÜV Rheinland	
CE Mark of a medical device	
Date of Manufacture	
Use by	
Lot Code	
Device Number	
Manufacturer	

Consult Instruction of Use	 https://www.bebdental.it/pro/en/instructions/
Do Not Re-use	
Caution	
Sterilised using irradiation Single sterile barrier system with protective packaging outside	
Do not use if package is damaged or open	
Unique Device Identifier	
Do not resterilize	
Non-Sterile Device	
Medical device	
Keep away from the sunlight	
Keep dry	

Table 2

Patient name or patient ID	
----------------------------	---

Date of implantation	
Name and Address of the implanting healthcare institution/provider	
Information website for patient	
Medical device	

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