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**DOLL MEMBRANE**

**General Information & Intended Use**

Non-reabsorbable titanium membranes are one of the most effective solutions for the stabilisation and protection of the graft material in the bone defect, as autologous bone and derivatives or artificial materials. Indeed, all types of titanium membranes are suitable for the treatment of all cases of guided bone regeneration of upper or lower jaw, especially in anticipation of subsequent rehabilitation with intraosseous implants.

In particular, titanium membranes and their screws can be used in the following applications:

- Increment or increase of bone tissue in the area of implant reception.
- Regeneration of residual bone defects due to extractions, cystectomies, chronic infectious processes.
- Increase of bone tissue, both in the vertical and labiobuccal direction, of the alveolar edge in order to obtain sufficient volume for the insertion of submerged implants.
- Improvement of healing conditions in implant areas.
- Dehiscence or fenestration of cortical plate with partial exposure of implants.
- Defects due to discrepancy between a post-extraction socket and the implant placed in it.

Titanium membranes have a microperforated surface, made from pure grade 2 titanium for medical use. They can be easily contoured, and cut using scissors, if necessary, to perfectly adapt them to the defect area so that the final margins are in contact with the residual healthy bone ridge on which they are stabilised.

Locking and osteosynthesis screws are made from titanium alloy.

**T-BARRIER Grids** should be used especially for large bone defect filling. They shall be fixed to the bone with osteosynthesis screws to keep the regeneration materials in place in the site.

**T-BARRIER grid with central hole, SIM and DOLL MEMBRANE** are suitable for the treatment of cases of guided bone regeneration with subsequent rehabilitation with intraosseous implants, preventing their migration in the maxillary sinus, thanks to the use of the dedicated locking screw.

**T-BARRIER grid with central hole** is more suitable in case of large bone defect filling.

**SIM** are suitable for maxillary sinus elevation.

**DOLL membranes** are specific for lateral bone defect filling.

**SUPPLY**

**Titanium membranes, locking and osteosynthesis screws are available in sterile primary packaging.**

**The primary packaging of titanium membranes is a blister containing one membrane or a locking screw.**

**Blisters are inserted into a protective cardboard box (secondary packaging).**

**Each box contains two blisters, traceability labels, package leaflet and implant card.**

**Traceability labels report the product code, product description, batch number and the relative UDI.**

**Codes available and related characteristics**

Line	Ref.	Dimensions			Cardboard box content	
		Width (mm)	Length (mm)	Thickness (mm)	Membranes	Locking screw
Grid	00532	23	28.5	0.13	2	No
Grid with central hole	00532/12	23	28.5	0.13	1	Yes
SIM	00532/3	11.8	19.5	0.13	1	Yes
	00532/4	13.8	20.1	0.13	1	Yes
	00532/5	17.6	19.7	0.13	1	Yes
DOLL	00532/1	8	18	0.13	1	Yes
	00532/2	8	24	0.13	1	Yes
	00532/10	7	10	0.13	1	Yes

## Warning

- Under law, only surgeons and certified dentists are allowed to use the medical devices in surgery. Any use by other persons is strictly prohibited. Therefore, titanium membranes and related accessories must be handled and placed only by qualified and specially trained personnel.
- B. & B. Dental assumes no responsibility for any loss or damage/injury to the patient or any serious health hazards caused through improper use, or any use by unqualified or inadequately trained personnel, that is not in accordance with this instruction manual.
- B. & B. Dental is not liable for any loss or damage/injury to the patient when accessories not manufactured by B. & B. Dental are used.
- Medical devices covered by these IFU must be placed in a hospital or clinical environment.
- Any prosthetic load on the regeneration area is strictly prohibited until it is healed.
- Please check "Table 1" for a better understanding of the symbols on medical device labels.
- Potential adverse events associated with the use of this product may include, but are not limited to, the following.
  - Dehiscence of the titanium membrane
  - Local or systemic infection

Proper selection of patient candidates, implant placement and post-operative care, combined with correct product use may help to mitigate these risks.

## General contraindications

Titanium membranes and related accessories shall not be placed in case of **non-compliant patient** and if there is any **intraoral contraindication**, or if there is clinical or radiographic evidence of local pathology. Chronic or acute infectious diseases, chronic maxillary sinus osteitis, vascular problems, systemic diseases shall also be evaluated.

## Absolute contraindications

**Known allergies to titanium, recent myocardial infarction, valvular prosthesis, severe kidney disorder, severe liver disorders, osteomalacia, insulin-resistant diabetes, recent high-dosage radiotherapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.**

## Relative contraindications

**Chemotherapy, moderate kidney disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged corticosteroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.**

## Precautionary information

- **Since the surgical techniques are highly specialised and complex procedures, it** is strongly recommended that practitioners shall complete a training course to learn established techniques of placement in oral implantology. Failure to seek proper training or follow instructions for use (IFU) may impact clinical success of implant, such as **failure of the implant, loss of bone, and postoperative complications, that could lead to pain for the patient and seriously endanger his or her health.**
- Surgical placement and consequent bone regeneration require complex dental procedures that may involve risks. It is the responsibility of the practitioner to inform the patient of any risk(s) and alternative treatments (including a "no treatment" option) prior to performing any procedure.
- Use of proper patient selection criteria is critical to clinical success:
  - Patients shall be adults, but, in some cases, they could be adolescents, over 13.
  - Ethnicity doesn't have any impact on the use of these medical devices.
  - The patient's current health status, medical and dental history have to be evaluated taking into account the contraindications, warnings and precautionary information.
  - Determine the location of all avoidable anatomical features prior to initiating any implant procedure, through X-rays, palpation, probing, and visual inspection of the proposed implant site.
  - Poor bone quality, poor patient oral hygiene habits, tobacco use, certain concomitant medications, and systemic disease (diabetes, etc.) may contribute to lack of osseointegration and subsequent implant failure.
- Titanium membrane, and locking and osteosynthesis screws are sterile and single-use medical devices. Do not sterilise and re-use. Any re-use and re-sterilisation of the medical devices are strictly prohibited as the functional performance and safety of the device expected by the manufacturer cannot be guaranteed, this would mean exposing the patient to the risk of loss of device performance, injury, pain and/or serious risks of contamination and cross-infections.

Proper use of surgical instrumentation is critical to clinical success:

- Any use and re-use of surgical instruments without an accurate and in-depth reading of their specific leaflets does not guarantee the functional performance and safety of the devices provided by the manufacturer, thus

exposing the patient to the risks of non-mechanical performance of the product and / or serious risks of contamination.

- All surgical instruments must be checked before each use regarding their state of cleanliness and sterilisation to ensure that they are always in operating condition. Do not use devices that show anomalies as they could compromise the safety and health of the patient.
- Inspect all surgical instruments before use. Using dull or worn instruments may cause the failure to correctly position the medical devices and consequent pain for the patient.
- Ensure there is a solid connection between any driver and handpiece and screw.
- Take adequate precautions to avoid that the patient inhales and/or ingests small components.
- Care must be taken when using instrumentation near other teeth to avoid damaging tooth structure.
- The pure titanium and titanium alloy used to manufacture titanium membranes and related accessories are not magnetic and thus should be acceptable for MRI procedures. However, it is important that patients advise medical professionals about their presence prior to imaging to ensure compatibility with the MRI equipment and to address potential artifacts, depending on the area of imaging.
- Instruct the patient to follow a routine post-surgery regimen as reported on the implant passport, according to the paragraph "**Post-surgery course**".
- Complications
  - Unsuccessful osseointegration is evidenced by possible infection, mobility or dehiscence.
  - Remove any exposed or not fixed membrane, and any **surrounding** visible granulation tissue, from the implant site as soon as possible, to prevent potential local complications.

### **Storage and use**

- All devices must be stored carefully in a clean and dry environment to keep the package intact.
- Devices should be stored at room temperature.
- Refer to individual product labels for special storage and handling conditions.
- Sterile products should be used before the expiration date printed on the label.
- Devices must not be used whenever the packaging is not intact.
- Contents are considered sterile unless the package is opened or damaged. Do not use products labelled "STERILE" if packaging has been damaged or opened prior to intended use.

### **Disposal**

Dispose of in accordance with the laws in force as special hospital waste. In particular, if the devices are contaminated with blood or fluids, it is necessary to use the appropriate containers and gloves to avoid any direct contact.

### **Surgical protocol for membrane placement**

#### ***Patient selection***

Careful case planning and appropriate patient selection criteria are critical and extremely important for successful clinical procedures. This should include assessment and review of the patient's current health status, medical and dental history, and the state and ability to maintain proper oral hygiene. Anatomical acceptability has to be evaluated, which implies adequacy of bone shall be determined by X-rays, palpation, probing, and thorough visual inspection of the proposed implant site. It is always recommended to perform a cone-beam CT scan to have precise anatomical references of thickness, height, density of bone and tissues and to define the defect. It is also essential to determine the location of all the anatomic structures that must be avoided before beginning any type of surgical procedure. Based on the assessment made, the type of regeneration to be performed shall be defined.

#### ***Implant selection***

The choice of the titanium membrane is important for long-term success. Evaluate the quality and type of bone in the regeneration site to determine the devices to use. In case of cone-beam CT scan, all the necessary anatomical references are immediately available for the evaluation.

#### ***Patient preparation***

Prepare the patient, following standard aseptic guidelines. Administer sufficient anaesthetic by local infiltration. In the mandible, inferior alveolar nerve block anaesthesia is not recommended in order to allow the patient to retain sensation and give feedback if necessary.

#### ***Surgical procedure***

As with any other surgical action, it is important for the procedure to be performed in an environment which is as sterile as possible. The surgical instrumentation is supplied sterile. It shall be disinfected, cleaned and sterilised before each use according to the specific IFU.

#### **Preparation of the site**

Make a full thickness mesiodistal incision along the surface of the crest. Using a periosteal elevator, leave the alveolar bone surface exposed in order to have adequate view and adequate working space for the surgical procedure.

The implant site must be prepared in a visible, clear field, so that the operator can prepare it with precision in every phase, thus allowing for proper positioning of the devices.

#### **Treatment and sterilisation**

Titanium membranes, and locking and osteosynthesis screws are supplied in sterile packages, and must be placed with cleaned and sterile surgical instruments. They have to be handled with titanium instruments and/or sterile, powder-free gloves to avoid contaminating the product or interfering with its biocompatibility.

#### **Membrane placement**

7. **According to the clinical case, the bone grafting material shall be positioned on the crest for filling in the defect. Meanwhile, the titanium membrane shall be contoured and possibly cut in order to perfectly fit onto the prepared area.**
8. Osteosynthesis screws have great penetration ability and preliminary site perforation is not generally required. It is, however, possible to perforate the site using the corresponding diameter size drill. **The perforation of the bone cortex has to be performed** by pumping the drill slightly up and down at a speed of **800-1000 RPM, under an abundant external irrigation with saline or, preferably, sterile distilled water. Furthermore, drilling must be intermittent both to avoid bone to heat up and to create a pumping effect that will help effective removal of bone tissue.**
9. **Then, the titanium membrane shall be placed on the prepared area and stabilised laterally with the use of the osteosynthesis screws, screwed with the dedicated spanner for contra-angle provided with cross-connection. In case of guided bone regeneration upon one-stage implant placement, the titanium membrane is fixed to the implant with a locking screw to prevent its movement. The locking screw, with a thread pitch of 1.8 mm, is suitable for DURA-VIT dental implant system. Always tighten the screws by turning clockwise until they are stable.**
10. **An X-ray check of the position of the implant is recommended before closing the flap.**
11. **Like all grafts, the primary stability of the grid is an indispensable condition for success. It is also essential to perfectly suture the mucosal flaps, so as to completely cover the grid, which is an important aspect to prevent dehiscence of the membrane and consequent failure.**
12. **Membranes and screws shall be removed after the bone regeneration period. This period must be established by the practitioner.**

#### **WARNING: Screwing instructions**

Please comply with the driving forces indicated here below:

- **Maximum 25 Ncm for placing the locking and osteosynthesis screws with a spanner for contra-angle.**

#### **Post-surgery course**

**The patient must be taught to follow the post-surgery course.**

The patient should be informed on how to behave immediately after surgery and in the subsequent days, which includes **cold compresses**, no physical exertion and/or sauna for the first 24 hours, no alcohol, no nicotine, no tea or coffee in the first days. The patient must be cautioned against chewing in the area where the membrane was inserted.

**Moreover, s/he should be advised that poor oral hygiene, tobacco smoking and systemic and general diseases (diabetes, rheumatoid arthritis, etc.) can contribute to poor bone integration and a consequent failure of the implant.**

**A prescription must be made for taking antibiotic an hour before the operation, and for 7 days after the surgery. The stitches must be removed after one week.**

Later, the patient shall maintain an appropriate and continuous oral hygiene, have regular check-ups and professional cleaning, and contact their dentist if s/he experiences tender points, notices loss of retention or membrane dehiscence. Patients shall advise medical professionals about the presence of the implanted medical devices prior to MRI. In most cases, these medical devices will not interfere with MRI.

**Customer information**

No person is authorised to provide any information which deviates from the information provided in this instruction sheet.

For further information and technical documents, please consult our website [www.bebdental.it](http://www.bebdental.it) at the “download” area.

**The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED.**

**In case any serious or suspected accidents occur following the use of our medical devices, you shall immediately report them to the Company and to your National Competent Authority.**

**Table 1**

CE Mark with the number of the Notified Body TÜV Rheinland	
CE Mark of a medical device	
Date of Manufacture	
Use by	
Lot Code	
Device Number	
Manufacturer	
Consult Instruction of Use	 <a href="http://www.bebdental.it/pro/en/instructions/">http://www.bebdental.it/pro/en/instructions/</a>
Do Not Re-use	
Caution	
Sterilized using irradiation Single sterile barrier system with protective packaging outside	
Do not use if package is damaged or open	
Unique Device Identifier	
Do not resterilize	

Non-Sterile Device	
Medical device	
Keep away from the sunlight	
Keep dry	

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