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WARNING: •The products covered by this instruction manual must be handled and positioned only by qualified personnel in a hospital and clinical environment. All the staff must be duly trained. Please be aware that proper handling by the staff prevents from physical hazards such as from sharps, as well as possible risks of infection. When necessary, qualified personnel shall explain to the patient the hygiene rules to be adopted and plan a regular outpatient check-up programme. •B. & B. Dental accepts no liability for any damage to the patient or serious health hazards caused by any improper use, not in accordance with the instructions for use and use by unqualified or inadequately trained personnel. •B. & B. Dental is not liable for malfunctions or breakages of products that have been modified in any way from their origin. •The use of instruments not manufactured by B. & B. Dental, releases the manufacturer from any liability for any damage to the patient. •For the specific intended use of the various families, please refer to the DEVICE LIST below. •In order to avoid overheating of the bone, which could compromise the outcome of implant surgery, the cutting speed, torque and all indications in the leaflet of the specific dental implant to be placed, which can be downloaded from the following link, must be strictly observed: <https://www.bebdental.it/pro/en/instructions/>. •For further information on surgical procedures, please refer to the 'Product Catalogue' and the company website. •Check 'Table 1' to better understand the symbols on medical device labels.

SPECIFIC WARNINGS FOR DISPOSABLE AND/OR STERILE DEVICES: •Any reuse of Disposable Devices is strictly prohibited, as the functional performance and the safety of the Device intended by the manufacturer cannot be guaranteed, exposing the patient to risks of mechanical failure of the product and/or serious risks of contamination and cross-infection. •Sterile devices must be used before the expiry date printed on the label. •Any re-sterilisation of Sterile Devices is strictly prohibited, as the functional performance and the safety of the Device as intended by the manufacturer cannot be guaranteed, exposing the patient to serious risk of contamination.

SPECIFIC WARNINGS FOR SINGLE-USE AND NON-STERILE DEVICES: Wash and sterilise disposable and non-sterile devices before use in accordance with the user's own washing and sterilisation instructions.

SPECIFIC WARNINGS FOR REUSABLE DEVICES: •All instruments must be washed, disinfected and sterilised the first time they are used and each time they are re-used in accordance with the following section '**Instructions for cleaning, disinfecting, drying and sterilising surgical instruments**'. •Any use and re-use of non-sterile surgical instruments without careful and thorough reading of this package leaflet and the catalogues below does not guarantee the functional performance and safety of the devices as intended by the manufacturer, exposing the patient to the risk of serious contamination and/or failure of the product's mechanical performance. •B. & B. Dental does not replace or repair instruments that have lost their functionality due to deterioration through normal use.

CONTRAINDICATIONS

They are strictly related to the specific indications reported on the IFU of dental implants and titanium membranes.

General contraindications

Do not use in case of non-compliant patient and if there is any intraoral contraindication, or if there is clinical or radiographic evidence of local pathology. Chronic or acute infectious diseases, chronic maxillary sinus osteitis, vascular problems, systemic diseases must also be evaluated. Check the quality and/or quantity of the bone tissue. Limited blood supply and previous infections that could delay or compromise healing.

Absolute contraindications

Supposed or confirmed allergy or sensitivity to steel or titanium, recent myocardial infarction, valvular prosthesis, severe kidney disorder, severe liver disorders, osteomalacia, insulin-resistant diabetes, recent high-dosage radiotherapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.

Relative contraindications

Chemotherapy, moderate kidney disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged corticosteroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Precautionary information

- Use of proper patient selection criteria is critical to clinical success:
 - Patients shall be adults, from the age of 18 years old.
 - Ethnicity doesn't have any impact on the use of these medical devices.
 - Due to the potential effect of the anaesthesia, do not use on pregnant women.
 - The patient's current health status, medical and dental history have to be evaluated taking into account the contraindications, warnings and precautionary information.
 - Determine the location of all the anatomical features to be avoided prior to initiating any implant procedure.
 - Patients with hyperfunction (bruxism) may be at increased risk for implant failure or fracture.
 - Adequacy of bone should be determined by X-rays, palpation, probing, and thorough visual inspection of the proposed implant site.

- Poor bone quality, poor patient oral hygiene habits, tobacco use, certain concomitant medications, and systemic disease (diabetes, etc.) may contribute to lack of osseointegration and subsequent implant failure.
- Proper use of surgical instrumentation is critical to clinical success:
 - The surgical instrumentation has been specifically designed for the dental implant lines by B. & B. Dental.
 - Any use and re-use of surgical instruments without an accurate and in-depth reading of their specific leaflets does not guarantee the functional performance and safety of the devices provided by the manufacturer, thus exposing the patient to the risks of non-mechanical performance of the product and / or serious risks of contamination and cross-infection.
 - All surgical instruments must be checked before each use regarding their state of cleanliness and sterilisation to ensure that they are always in operating condition. Do not use devices that show anomalies as they could compromise the safety and health of the patient.
 - Inspect all surgical instruments before use. Using dull or worn instruments may cause implant failure.
 - Ensure there is a solid connection between any driver or drill, handpiece or any eventual accessory.
 - Take adequate precautions to avoid that the patient inhales and/or ingests small components.
 - Care must be taken when using instrumentation near other teeth to avoid damaging tooth structure.
- The accessories for dental implants and titanium membranes are manufactured from not magnetic materials and thus should be acceptable for MRI procedures. However, it is important that patients advise medical professionals about the presence of dental implants prior to imaging to ensure compatibility with the MRI equipment and to address potential artifacts, depending on the area of imaging.

Storage and use

- All devices must be stored carefully in a clean and dry environment to keep the package intact.
- Devices should be stored at room temperature.
- Refer to individual product labels for special storage and handling conditions.
- Sterile products must be used before the expiration date printed on the label.
- Devices must not be used whenever the packaging is not intact.
- Contents are considered sterile unless the package is opened or damaged. Do not use products labelled "STERILE" if packaging has been damaged or opened prior to intended use.

Disposal

Dispose of in accordance with the laws in force as special hospital waste. In particular, if the devices are contaminated with blood or fluids, it is necessary to use the appropriate containers and gloves to avoid any direct contact.

ACCESSORIES FOR DENTAL IMPLANTS

General indications: Accessories for dental implants have been specially designed for the B. & B. Dental implant lines. The user has the duty to determine which device is suitable and appropriate for the patient and the specific clinical case. For the specific intended use of the various families, please refer to the DEVICE LIST below.

DEVICE LIST

A) COVERING SCREWS - Intended use: Covering the implant cavity during osseointegration - Description: The threaded part must be screwed into the internal thread of the implant, while the head closes the implant cavity. Manufactured from Grade 5 Titanium. Available in different sizes and heights; Disposable and sterile devices. Present in implant packaging or supplied in individual packages.

B) HEALING SCREWS - Intended use: Soft tissue healing in order to create an ideal fit for prosthetic components. - Description: They fit into the implant thanks to the threaded part. Manufactured from Grade 5 Titanium and peek and available in different sizes and heights to adapt to various gingival thickness. They are disposable and sterile devices, supplied in individual packages, except for the transmucosal peek screw which is included in the packaging of some implants.

C) TRANSFER - Intended use: Accessories for copying the implant position onto the laboratory model. - Description: Products made of Grade 5 Titanium and Delrin. Disposable and sterile devices. Supplied in single packs. - Family detail: for open tray, for closed tray, pull-off and for digital intraoral scans.

D) ABUTMENTS- Intended use: Prosthetic options to produce superstructures and dental prostheses to ensure optimum resistance to masticatory load. - Description: Produced in titanium, peek, Plexiglas, cobalt chrome, medical polymers. Available in different shapes and/or transmucosal heights to better suit the aesthetic needs of the patient. Disposable devices, supplied in single, non-sterile packages, the only exception being spherical abutments, which are sold sterile. Connection: Tapered fit with the implant and fastening by means of the prosthetic screw included in the package. Only the spherical abutments fit into the implant by means of the threaded part. Tightening must be carried out according to the instructions in the section "*Screwing and tightening instructions*". - Family detail: •Temporary abutments •UCLA abutments •Straight and angled titanium abutments •Digital components •Multi use abutments (MUA) •Conical abutments •Flat abutments •Spherical anchoring system.

E) PROSTHETIC SCREWS - Description: Manufactured from grade 5 titanium. Available in different models and heights. Disposable devices supplied in individual packages or together with their own component. - Family detail: •Prosthetic screws, supplied non-sterile • Transfer screws, sold sterile - Intended use Tightening of prosthetic components and transfers. - **ATTENTION:** *Observe the tightening instructions given in the section 'Screwing and tightening instructions'.* Please refer to the 'Product Catalogue' for further screwing instructions.

F) ANALOGUES - Intended use: Defining the position of the device in the laboratory model. - Description: Manufactured from grade 5 titanium. Available in different models. Disposable devices supplied in single packs.

G) CAPS - Intended use: Accessories for the connection of mobile prostheses made with a spherical anchoring system. - Description: Grade 5 titanium with internal medical polymer cap or O-ring. Available in different models. Disposable devices. Supplied in single packs.

H) IMPRESSION COPING FOR MINI IMPLANTS - - Intended use: Accessories for copying the position of the mini-implant onto the laboratory model. - Description: Products made of Derlin, PEEK and titanium; Disposable, non-sterile devices; Supplied in single packs.

I) SCAN CAP - Intended use: Accessories for intraoral digital scanning. - Description: Manufactured from Peek and supplied with grade 5 Titanium connection screws. Available in different models. Disposable devices. Supplied in single packs.

J) KEYS - Intended use: Screwing of medical devices. - Description: Manufactured from medical grade stainless steel; Reusable instruments; Available in different diameters, lengths and models; Supplied in individual, sealed and NON-STERILE packages. - Family detail: **1) for contra-angle:** •Implant insertion keys •Prosthetic screwdriver •Connectors **2) for ratchet:** •Implant Insertion keys •Prosthetic Screwdrivers •Spherical abutment drivers •Extractor keys •Mounters •Mounter drivers **3) manual:** •Manual keys •Prosthetic screwdrivers •Butterfly keys for mini-implants •Spherical abutment drivers •Positioners for abutments •Extractors • Ratchet key adaptors **4) ratchets:** •Torque ratchets •Universal ratchets **5) Prosthetic kits.**

Screwing and tightening instructions

- 1) **For the insertion of all implants with contra-angle keys**, 35 Ncm must not be exceeded, otherwise, **with ratchet keys**, the tightening torque values which must be strictly adhered to are as follows: For implants with a diameter of 3.5 mm or more do not exceed 50 Ncm. For SLIM implants do not exceed 35 Ncm.
- 2) 45 Ncm must not be exceeded when screwing in mini-implants
- 3) The screwing force of the **cover screws** must not exceed 10 Ncm;
- 4) The screwing force of the **healing and transmucosal screws** must not exceed 20 Ncm.
- 5) The following tightening torques must be strictly adhered to in order to prevent **prosthetic screw** from breaking: 25 Ncm, if the prosthesis is to be clamped on implants with a diameter of 3.5 mm or more. 20 Ncm, if the prosthesis is to be clamped on SLIM implants.
- 6) For screwing **abutments (multi-use abutments, flat and spherical abutments)**, the following screwing torques must be strictly observed: 25 Ncm, if abutments are to be fixed on implants with a diameter of 3.5 mm or more. 20 Ncm, to fix abutments on SLIM implants and abutments for FLAT for implants with a diameter of 3.5 mm or more. 15 Ncm, if abutments for MUA or FLAT are to be fixed for SLIM implants.
- 7) **Manual** - Apply a force appropriate to the clinical case.
- 8) **Ratchets** - **WARNING:** Instructions for use: When screwing (IN side), the ratchets must be used in a horizontal position; therefore, the force applied will be lateral and not apical. During unscrewing (OUT side), the ratchets perform their function even when they are not aligned.

Fixed ratchet: The fixed ratchet, unlike the torque ratchet, allows implants to be screwed or unscrewed without adjusting the force applied. It should only be used by experienced and qualified medical personnel. •Replacing the ratchet wheel: Pull back the pin on both sides with your thumb and forefinger, then pull out or insert the wheel, as required. Then release the pin. Please note: the wheel can only be inserted from one side. •Correct operation: press only at the end of the handle and only with one finger.

Peek torque ratchet: The Peek ratchet allows implants to be screwed or unscrewed with direct tightening adjustment. The clamping force is indicated on the metal tab on the inside of the handle. The latter moves according to the force used by the doctor, showing its value.

Steel torque ratchet: Setting the prosthesis-tightening function: The tightening can be continuously adjusted with the adjusting nut via the spring. The setting can be read on the scale of the scale capsule. Setting the lock-in function: Turn the adjusting nut to graduation. Don't screw too tight! Store the ratchet with the spring unloaded! Replacing the ratchet wheel: Pull back the pin on both sides with your thumb and forefinger, then pull out or insert the wheel, as required. Then release the pin. Please note: the wheel can only be inserted from one side.

Correct functioning: Tightening can only be adjusted by turning the head of the adjusting nut. When the desired torque is reached, the scale capsule bends around the axis in the ratchet head. The release is audible and perceptible. *Stop pressing when torque is reached - the ratchet could be damaged.* Loosening the adjusting nut returns the ratchet to its initial position.

Reuse and disassembly instructions

Fixed ratchet: After use, remove the ratchet wheel by pushing it outwards and unscrew the key head so that it separates from the handle. This operation does not require any tools. Proceed with cleaning and disinfection as described in the next section, then reassemble the ratchet before sterilisation.

Peek torque ratchet: After use, remove the ratchet wheel by pushing it outwards. This operation does not require any tools. Proceed with cleaning and disinfection as described in the next section, then reassemble the ratchet before sterilisation.

Steel torque ratchet: After use, disassemble the ratchet into its parts as shown in the image below. This operation does not require any tools. Proceed with cleaning and disinfection as described in the next section, then reassemble the ratchet before sterilisation.

K) SURGICAL KIT - Intended Use: Organisation and housing of instruments to perform surgical procedures. - Description: The box is made of RADEL with silicone inserts. The size, composition and content vary depending on the implant line and the surgical technique to be adopted according to the clinical case. Supplied in cases, containing a single, non-sterile surgical kit. Reusable devices. For further information, please refer to the 'Product Catalogue' and the company website.

ACCESSORIES FOR TITANIUM MEMBRANES

General indications: The user has the duty to determine which device is suitable and appropriate for the patient and the specific clinical case. For the specific intended use of the various families, please refer to the DEVICE LIST below.

DEVICE LIST

A) OSTEOSYNTHESIS SCREWS - Intended use: Fixation and stabilisation of titanium membranes and custom-made implants. - Description: They have a threaded body, with a wide or narrow loop, which must be screwed into the bone. Manufactured from Grade 5 Titanium. Available in different sizes and heights. Disposable and sterile devices. Supplied in single packs. - Family detail: •Osteosynthesis screws •Fixing pins.

B) LOCKING SCREWS - Intended use: Fixing and stabilisation of titanium grid membranes with central hole, DOLL and SIM. - Description: The threaded body must be screwed into the implant and is specially designed for titanium membranes and B. & B. Dental implants.; Manufactured from Grade 5 Titanium; Available in different sizes and heights; Disposable and sterile devices; Supplied in single packs.

C) KEYS - Intended use: Screwing of medical devices. - Description: Manufactured from medical grade stainless steel; Reusable instruments; Available in different diameters, lengths and models; Supplied in individual, sealed and NON-STERILE packages. - Family detail: **1) for contra-angle:** •Keys for osteosynthesis screws •Connectors **2) manual:** •Manual keys for osteosynthesis screws.

INSTRUCTIONS FOR CLEANING, DISINFECTING, DRYING AND STERILISING SURGICAL INSTRUMENTS

It is recommended to carefully follow the instructions below to achieve proper and effective cleaning and sterilisation of the instrumentation.

WARNINGS: •Strictly adhere to the steps described below and carry them out sequentially, without skipping any steps. •Always use personal protective equipment (such as gloves, mask, goggles, etc.) at all stages. •B. & B. Dental has validated its surgical instrumentation for up to 30 reuses. •The operator must check the functionality of the equipment after each use. •The devices must not be placed in contact with other objects during washing. •Wash any surgical box and instruments separately. •Do not use metal-bristled brushes or steel wool for cleansing. •Do not use chlorine-based disinfectants under any circumstances. •Use approved and validated washer disinfectors according to ISO 15883. •Instruments should only be reinserted, after possible reassembly, into the surgical box only after all cleaning steps have been completed. •Chemical sterilisation is not recommended as this procedure may damage the plastic surface of the surgical boxes and may oxidise the devices. •Do not use dry sterilisation as the high temperature (approx. 180 °C) may melt the plastic surgical cassette. •Do not sterilise in glass bead sterilisers (temperature above 300 °C). •Subject the autoclave, ultrasonic bath and washer disinfectant to regular maintenance and checks as stipulated in the user manual provided by the manufacturer. •It is recommended to validate your autoclave at full and worst load according to ISO 17665-1.

CLEANING STEPS - MANUAL OPTION

1. PRE-TREATMENT: Dismantle multiple instruments into their respective parts (where applicable). Remove the coarsest impurities from the instruments by rinsing under running water for at least two minutes. This should be done within a maximum of 2 hours after use.

2. CLEANSING: •Immerse the devices immediately after rinsing under running water (pre-treatment) in a basin of detergent solution for at least one hour to aid cleansing and prevent any residual blood from coagulating. •The devices must be completely immersed in the solution and the bowl must be shaken several times to facilitate the process. •We recommend using a peracetic acid solution, such as 'SEKUSEPTTM PULVER CLASSIC' in a 2% concentration. •It is particularly necessary to clean well inside the cavities of the devices. •Wipe any residual blood, saliva, tissue or bone from the instruments immediately after surgery with a soft-bristled toothbrush under running water for at least 2 minutes. Residues adhering to instruments can cause corrosion as they dry out.

3. DISINFECTION: •For thorough cleaning, place the devices in the ultrasonic machine and leave them immersed for a minimum of 30 minutes to a maximum of one hour. To avoid corrosion of surgical instruments, this liquid must be suitable for disinfecting metal and plastic instruments. We recommend a disinfectant based on peracetic acid, such as 'SEKUSEPTTM PULVER CLASSIC' at 4%. • Finally, carry out a thorough and abundant rinsing preferably with sterile water and subsequent drying as described in the next section. At the end of the cleaning procedure, the instruments must be free of organic residues and other impurities.

4. DRYING PHASE

Dry the cleaned devices in a hot air oven at 60 °C for 2 hours. It is recommended to carefully follow the operating instructions provided by the supplier of such a dryer. Carefully check that the devices are completely dry, as any vapour residue can trigger corrosive processes that will damage the instrument over time. Then proceed to the sterilisation phase.

AUTOMATED CLEANING, DISINFECTION AND DRYING OPTION

•To carry out automated cleaning, disinfection and drying, follow the instructions provided by the manufacturer of the washer disinfectant in use. •B. & B. Dental validated the cycle with a washer disinfectant model AT-OS AWD655-10 (DIN 10) that is programmable and equipped with data integrity validated SW. The process is described below.

•Place the devices and any surgical box in the appropriate basket suitable for mechanical washing, close it and load it into the washer disinfectant. •It is recommended to avoid contact of the devices with each other during power washing. •Start the appropriate washing, disinfection and drying programme for the devices. •It requires a pre-rinse with demineralised water for 3 minutes. •Cleaning should be carried out for 5 minutes at a temperature of 55°C, followed by rinsing for 2 minutes at a temperature of 40°C. •We recommend using the concentrated liquid alkaline detergent Neodisher MediClean Forte (0.5%) and the alkaline residue neutraliser Neodisher Z (0.1%). •Final rinse with demineralised water for 2 minutes. •The cycle continues with automated disinfection at 90°C for 5 min. •Take the devices from the washer disinfectant at the end of the drying programme and proceed to the sterilisation phase.

PHASES OF STERILISATION

1.PACKAGING BEFORE STERILISATION: •Disinfected, clean and dry devices must be placed neatly in the surgical cassette, sealed in special sterilisable pouches validated in accordance with ISO 11607-1. Alternatively, they can be placed individually in paper (60 g/m²) and plastic film (PET/PP) envelopes. •Paper envelopes (60 g/m²) and plastic film (PET/PP) must be declared steam sterilisable by the supplier. •Paper and plastic film envelopes must be sealed on all four sides. •B. & B. Dental suggests using the same sealer model and sealing parameters used during validation: sealer model LUX SEAL M9025/4; sealing time: 5 seconds; sealing temperature: 150 °C.

2.STERILISATION: •Sterilisation must be carried out in a steam autoclave at a temperature of 121 °C for a minimum of 15 minutes. • After sterilisation, leave to dry for 20 minutes. •Store the box or individually wrapped devices in a dry environment and within their own envelope for up to 6 months from the sterilisation date. •B. & B. Dental suggests using an autoclave model equivalent to the one used in the validation phase: model "SERENA 18" (Manufacturer Reverberi Srl).

For further information and technical documents, please consult our website www.bebdental.it, especially the "download" page.

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED, through the following website: <https://ec.europa.eu/tools/eudamed>. The SSCP is linked to the following BASIC UDI-DI:

- **805228249BIPHASICIMPLAN6P (biphasic implants);**
- **805228249MONOPHASICIMPLANB6 (monophasic implants);**
- **805228249MINIIMPLANTS DP (mini implants);**
- **805228249TEMPORARYABTCL (titanium temporary abutments);**
- **805228249PTEMPORARYABT7G (castable and peek temporary);**
- **805228249FINALABUTMENTXU (final abutments);**
- **805228249SPHERICALABTBZ (spherical abutments);**
- **805228249HEALINGDEVICESWA (healing screws, cover and transmucosal screws);**
- **805228249OSTEOSCREWPINAM (osteosynthesis and locking screws and pins).**

In the event that serious or suspected accidents occur following the use of our medical devices, you shall report them immediately to the Company and to your National Competent Authority.

Table 1

CE Mark with the number of the Notified Body TÜV Rheinland	
CE Marking of a Medical Device	
Date of Manufacture	
Use by	
Lot Code	
Device Number	
Manufacturer	
Consult Instruction of Use	 https://www.bebdental.it/pro/en/instructions
Do Not Re-use	
Attenzione / Caution	
Sterilized using irradiation Single sterile barrier system	
Do not use if package is damaged or open	
Unique Device Identifier	
Do not resterilize	

Non-Sterile Device	
Medical device	
Keep away from the sunlight	
Keep dry	

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