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WARNING: •The products covered by this instruction manual must be handled and positioned only by qualified personnel in a hospital and clinical environment. All the staff must be duly trained. Please be aware that proper handling by the staff prevents from physical hazards such as from sharps, as well as possible risks of infection. When necessary, qualified personnel shall explain to the patient the hygiene rules to be adopted and plan a regular outpatient check-up programme. •B. & B. Dental accepts no liability for any damage to the patient or serious health hazards caused by any improper use, not in accordance with the instructions for use and use by unqualified or inadequately trained personnel. •B. & B. Dental is not liable for malfunctions or breakages of products that have been modified in any way from their origin. •The use of instruments not manufactured by B. & B. Dental, releases the manufacturer from any liability for any damage to the patient. •For the specific intended use of the various families, please refer to the DEVICE LIST below. •In order to avoid overheating of the bone, which could compromise the outcome of implant surgery, the cutting speed, torque and all indications in the leaflet of the specific dental implant to be placed, which can be downloaded from the following link, must be strictly observed: <https://www.bebdental.it/pro/en/instructions/>. •For further information on surgical procedures, please refer to the 'Product Catalogue' and the company website. •Check 'Table 1' to better understand the symbols on medical device labels.

SPECIFIC WARNINGS FOR REUSABLE DEVICES: •All instruments must be washed, disinfected and sterilised the first time they are used and each time they are re-used in accordance with the following section 'Instructions for cleaning, disinfecting, drying and sterilising surgical instruments'. •Any use and re-use of non-sterile surgical instruments without careful and thorough reading of this package leaflet and the catalogues below does not guarantee the functional performance and safety of the devices as intended by the manufacturer, exposing the patient to the risk of serious contamination and/or failure of the product's mechanical performance. •B. & B. Dental does not replace or repair instruments that have lost their functionality due to deterioration through normal use..

CONTRAINDICATIONS

They are strictly related to the specific indications reported on the IFU of dental implants.

General contraindications

Do not use in case of non-compliant patient and if there is any intraoral contraindication, or if there is clinical or radiographic evidence of local pathology. Chronic or acute infectious diseases, chronic maxillary sinus osteitis, vascular problems, systemic diseases must also be evaluated. Check the quality and/or quantity of the bone tissue. Limited blood supply and previous infections that could delay or compromise healing.

Absolute contraindications

Supposed or confirmed allergy or sensitivity to steel or titanium, recent myocardial infarction, valvular prosthesis, severe kidney disorder, severe liver disorders, osteomalacia, insulin-resistant diabetes, recent high-dosage radiotherapy, severe chronic alcoholism, drug abuse, illnesses and tumours in the terminal phase.

Relative contraindications

Chemotherapy, moderate kidney disorders, moderate liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding or motivation, AIDS, HIV-seropositivity, prolonged corticosteroid use, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Precautionary information

Use of proper patient selection criteria is critical to clinical success:

- Patients shall be adults, at least aged 18 years old.
- Ethnicity doesn't have any impact on the use of these medical devices.
- Due to the potential effect of the anaesthesia, do not use on pregnant women.
- The patient's current health status, medical and dental history have to be evaluated taking into account the contraindications, warnings and precautionary information.
- Determine the location of all avoidable anatomical features prior to initiating any implant procedure, through X-rays, palpation, probing, and visual inspection of the proposed implant site.
- Adequacy of bone should be determined by X-rays, palpation, probing, and thorough visual inspection of the proposed implant site.
- Poor bone quality, poor patient oral hygiene habits, tobacco use, certain concomitant medications, and systemic disease (diabetes, etc.) may contribute to lack of osseointegration and subsequent implant failure.

Proper use of surgical instrumentation is critical to clinical success:

- The surgical instrumentation has been specifically designed for the dental implant lines by B. & B. Dental.

- Any use and re-use of surgical instruments without an accurate and in-depth reading of their specific leaflets does not guarantee the functional performance and safety of the devices provided by the manufacturer, thus exposing the patient to the risks of non-mechanical performance of the product and / or serious risks of contamination and cross-infections.
- All surgical instruments must be checked before each use regarding their state of cleanliness and sterilisation to ensure that they are always in operating condition. Do not use devices that show anomalies as they could compromise the safety and health of the patient.
- Inspect all surgical instruments before use. Using dull or worn instruments may cause the failure of a correct positioning of the medical devices and consequent pain for the patient.
- Ensure there is a solid connection between any drill, handpiece and the eventual accessory.
- Take adequate precautions to avoid that the patient inhales and/or ingests small components.
- Care must be taken when using instrumentation near other teeth to avoid damaging tooth structure.

An efficient and atraumatic implant site preparation is created through a procedure relying on a gradual drilling technique. The whole stage of bone tissue drilling must be performed under an abundant external irrigation with saline solution or, preferably, sterile double-distilled water. Furthermore, drilling must be intermittent both to avoid bone to heat up and to create a pumping effect that will help effective removal of bone tissue.

Storage and use

- All devices must be stored carefully in a clean and dry environment to keep the package intact.
- Devices should be stored at room temperature.
- Refer to individual product labels for special storage and handling conditions.
- Devices must not be used whenever the packaging is not intact.

Disposal

Dispose of in accordance with the laws in force as special hospital waste. In particular, if the devices are contaminated with blood or fluids, it is necessary to use the appropriate containers and gloves to avoid any direct contact.

DEVICE LIST

A) DENTAL DRILLS - Intended use: Preparation of the surgical site for subsequent positioning of the implants. Dental drills penetrate the soft tissues, perforate the jaw bone to create the surgical site suitable to host the implant. - Description: Manufactured from medical grade stainless steel. Available in different diameters, lengths and designs. Supplied in individual, NON-STERILE packages. - Family detail: • Lance drills • Rounded drills • Countersink drills • Depth and conical drills • Wide drills • Guided surgery drills • Mini implant or preparation drills • Pterygo implant drills • Mucotomes • Bone mills • Pin drills • Osteosynthesis screw drills • Bone taps - WARNING: • Dental drills must be used in combination with the contra-angle. Manual operation is not possible. • The drill of diameter 1.1 mm is considered SINGLE USE and any reuse of Disposable Devices is strictly prohibited, as the functional performance and safety of the Device intended by the manufacturer cannot be guaranteed, exposing the patient to risks of mechanical failure of the product and/or serious risks of contamination. Before using it, the procedure described in the paragraph "Instructions for cleaning, disinfection, drying and sterilization of surgical instruments" must be carried out.

B) BONE COMPACTING TOOLS - Family detail: • Compactors-Expanders • Push pins • Push screws - Intended use: Bone expansion and compaction for surgical site preparation. - Description: Manufactured from medical grade stainless steel. Available in different diameters, lengths and designs. Supplied in individual, NON-STERILE packages.

C) ACCESSORIES FOR DRILLS - Intended use: Accessories to make surgical site preparation easier and safer - Description: Manufactured from medical grade stainless steel or titanium. Available in different diameters and lengths. Supplied in single packs and NON-STERILE. - Family detail: • Extenders • Metal stopper • Converters • Parallelometer • Depth gauges • Crestal pins • Lateral pins • Arch guide • Spacer • Spacer pin.

D) SURGICAL KIT - Intended Use: Organisation and housing of instruments to perform surgical procedures. - Description: The box is made of RADEL with silicone inserts. The size, composition and content vary depending on the implant line and the surgical technique to be adopted according to the clinical case. Supplied in cases, containing a single, non-sterile surgical kit. For further information, please refer to the 'Product Catalogue' and the company website.

INSTRUCTIONS FOR CLEANING, DISINFECTING, DRYING AND STERILISING SURGICAL INSTRUMENTS

It is recommended to carefully follow the instructions below to achieve proper and effective cleaning and sterilisation of the instrumentation.

WARNINGS: • Strictly adhere to the steps described below and carry them out sequentially, without skipping any steps. • Always use personal protective equipment (such as gloves, mask, goggles, etc.) at all stages. • B. & B. Dental has validated its surgical instrumentation for up to 30 reuses. • The operator must check the functionality of the

equipment after each use. •The devices must not be placed in contact with other objects during washing. •Wash any surgical box and instruments separately. •Do not use metal-bristled brushes or steel wool for cleansing. •Do not use chlorine-based disinfectants under any circumstances. •Use approved and validated washer disinfectors according to ISO 15883. •Instruments should be reinserted, after possible reassembly, into the surgical box only after all cleaning steps have been completed. •Chemical sterilisation is not recommended as this procedure may damage the plastic surface of the surgical boxes and may oxidise the devices. •Do not use dry sterilisation as the high temperature (approx. 180 °C) may melt the plastic surgical cassette. •Do not sterilise in glass bead sterilisers (temperature above 300 °C). •Subject the autoclave, ultrasonic bath and washer disinfectant to regular maintenance and checks as stipulated in the user manual provided by the manufacturer. •It is recommended to validate your autoclave at full and worst load according to ISO 17665-1.

CLEANING STEPS - MANUAL OPTION

1. PRE-TREATMENT: Dismantle multiple instruments into their respective parts (where applicable). Remove the coarsest impurities from the instruments by rinsing under running water for at least two minutes. This should be done within a maximum of 2 hours after use.

2. CLEANSING: •Immerse the devices immediately after rinsing under running water (pre-treatment) in a basin of detergent solution for at least one hour to aid cleansing and prevent any residual blood from coagulating. •The devices must be completely immersed in the solution and the bowl must be shaken several times to facilitate the process. •We recommend using a peracetic acid solution, such as 'SEKUSEPTTM PULVER CLASSIC' in a 2% concentration. •It is particularly necessary to clean well inside the cavities of the devices. •Wipe any residual blood, saliva, tissue or bone from the instruments immediately after surgery with a soft-bristled toothbrush under running water for at least 2 minutes. Residues adhering to instruments can cause corrosion as they dry out.

3. DISINFECTION: •For thorough cleaning, place the devices in the ultrasonic machine and leave them immersed for a minimum of 30 minutes to a maximum of one hour. To avoid corrosion of surgical instruments, this liquid must be suitable for disinfecting metal and plastic instruments. We recommend a disinfectant based on peracetic acid, such as 'SEKUSEPTTM PULVER CLASSIC' at 4%. • Finally, carry out a thorough and abundant rinsing preferably with sterile water and subsequent drying as described in the next section. At the end of the cleaning procedure, the instruments must be free of organic residues and other impurities.

4. DRYING PHASE

Dry the cleaned devices in a hot air oven at 60 °C for 2 hours. It is recommended to carefully follow the operating instructions provided by the supplier of such a dryer. Carefully check that the devices are completely dry, as any vapour residue can trigger corrosive processes that will damage the instrument over time. Then proceed to the sterilisation phase.

AUTOMATED CLEANING, DISINFECTION AND DRYING OPTION

•To carry out automated cleaning, disinfection and drying, follow the instructions provided by the manufacturer of the washer disinfectant in use. •B. & B. Dental validated the cycle with a washer disinfectant model AT-OS AWD655-10 (DIN 10) that is programmable and equipped with data integrity validated SW. The process is described below.

•Place the devices and any surgical box in the appropriate basket suitable for mechanical washing, close it and load it into the washer disinfectant. •It is recommended to avoid contact of the devices with each other during power washing. •Start the appropriate washing, disinfection and drying programme for the devices. •It requires a pre-rinse with demineralised water for 3 minutes. •Cleaning should be carried out for 5 minutes at a temperature of 55°C, followed by rinsing for 2 minutes at a temperature of 40°C. •We recommend using the concentrated liquid alkaline detergent Neodisher MediClean Forte (0.5%) and the alkaline residue neutraliser Neodisher Z (0.1%). •Final rinse with demineralised water for 2 minutes. •The cycle continues with automated disinfection at 90°C for 5 min. •Take the devices from the washer disinfectant at the end of the drying programme and proceed to the sterilisation phase.

PHASES OF STERILISATION

1. PACKAGING BEFORE STERILISATION: •Disinfected, clean and dry devices must be placed neatly in the surgical cassette, sealed in special sterilisable pouches validated in accordance with ISO 11607-1. Alternatively, they can be placed individually in paper (60 g/m²) and plastic film (PET/PP) envelopes. •Paper envelopes (60 g/m²) and plastic film (PET/PP) must be declared steam sterilisable by the supplier. •Paper and plastic film envelopes must be sealed on all four sides. •B. & B. Dental suggests using the same sealer model and sealing parameters used during validation: sealer model LUX SEAL M9025/4; sealing time: 5 seconds; sealing temperature: 150 °C.

2. STERILISATION: •Sterilisation must be carried out in a steam autoclave at a temperature of 121 °C for a minimum of 15 minutes. • After sterilisation, leave to dry for 20 minutes. •Store the box or individually wrapped devices in a dry environment and within their own envelope for up to 6 months from the sterilisation date. •B. & B. Dental suggests using an autoclave model equivalent to the one used in the validation phase: model "SERENA 18" (Manufacturer Reverberi Srl).

For further information and technical documentation, please consult the download section at www.bebdental.it.

In case any serious or suspected accidents occur following the use of our medical devices, you shall immediately report them to the Company and to your National Competent Authority.

Table 1: Symbols in use

CE Mark with the number of the Notified Body TÜV Rheinland	
CE Marking of a Medical Device	
Date of Manufacture	
Use by	
Lot Code	
Device Number	
Manufacturer	
Consult Instruction of Use	 https://www.bebdental.it/pro/en/instructions/
Do Not Re-use	
Caution	
Do not use if package is damaged or open	
Unique Device Identifier	
Non-Sterile Device	

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